


# BE1000 HW1 - Spec Lab Report


Maxx Yung

1. The stock solution was prepared at  $0.03\text{g}/\text{L}$ , while the molecular weight is  $332.31\text{g}/\text{mol}$ . That means  $\frac{0.03\text{g}}{332.31\text{g}/\text{mol}} = 0.0000902772\text{mol}$ . Using the formula to calculate molarity, which is  $M = \frac{\text{mol}}{\text{L}}$ , we get  $\frac{0.0000902772}{1} = 0.0000902772\text{ M}$ . After converting from  $M$  to  $\mu\text{M}$ , we get  $90.2772\mu\text{M}$ .
2. Given final diluted concentrations in  $\mu\text{M}$ :
  1. 25 -  $275\mu\text{L}$  of stock and  $725\mu\text{L}$  of water.
  2. 20 -  $220\mu\text{L}$  of stock and  $780\mu\text{L}$  of water.
  3. 15 -  $165\mu\text{L}$  of stock and  $835\mu\text{L}$  of water.
  4. 10 -  $110\mu\text{L}$  of stock and  $890\mu\text{L}$  of water.
  5. 5 -  $55\mu\text{L}$  of stock and  $945\mu\text{L}$  of water.
3. For the  $10\mu\text{M}$  sample, we needed  $110\mu\text{L}$  of stock and  $890\mu\text{L}$  of water. Then, for the  $1000\mu\text{L}$  pipette for water, the number was 089. For the  $200\mu\text{L}$  pipette for the stock solution, the number was 110.



- 4.
5. Max absorbance according to our data is: 485.0nm

### Absorbance (at 485.0 nm) vs. Concentration



6.

The extinction coefficient is the slope of the line because  $A = \epsilon lc$ , but  $l = 1$ , so  $\frac{A}{c} = \epsilon$ . That means the extinction coefficient is  $0.0287 \text{ mol}^{-1} \text{ cm}^{-1}$ .

7. Our unknown sample had an absorbance value of 0.463741054. Based on the standard curve formula,  $A = \epsilon lc \rightarrow 0.463741054 = 0.0287 \times 1 \times c$ , so  $c = 16.158224878 = 16.16 \mu\text{M}$